

ABT 3 STEP THAW KIT



ABT 3 STEP THAW KIT is made up of 3 different solutions used for rehydrating thawed embryos after they have been frozen in 10% Glycerol. Each solution contains Glycerol and Sucrose at different levels to more slowly thaw embryos and protect them from rapid changes in cell volume that might occur with faster thawing procedures.

ABT 3 Step Thaw Kits start with ultra-purified Type 1 water and includes Gentamycin as an antibiotic and certified New Zealand sourced BSA. All ABT products are manufactured in a cleanroom environment and go thru strict quality control practices, ensuring consistency and quality from lot to lot.

- **Thaw 1:** Contains 0.5 M Sucrose, 5% Glycerol
- **Thaw 2:** Contains 0.5 M Sucrose, 2.5% Glycerol
- **Thaw 3:** Contains 0.6M Sucrose, 0% Glycerol
- 100% traceable components including certified New Zealand sourced BSA
- Filtered to 0.1uM; no additional filtering required
- Tight quality control standards equal consistent results
- Available in durable EVA bags and 7ml LDPE Single Use Dropper Bottle
- Made in the USA with two decades of experience

PRODUCT #	SIZE
ABT048	50 ml Bag
ABT248	20 ml Bag
ABT748	7 ml Single Use Bottle

Storage: 2 - 8° C
Shelf Life: 18 months from date of manufacture

ABOUT ABT 360

Located in Pullman, Washington, USA, ABT 360 is owned and operated by Kathy Shearer who has 25 years of experience manufacturing embryo transfer media and Alyssa Hood who has assisted in and managed that same manufacturing facility since 1995. Formed with integrity and experience, ABT 360 has what it takes to support all of your embryo transfer media needs. We look forward to continuing our relationships with the veterinary reproduction community and we will strive to provide quality and consistency in our products, supported by exceptional customer service.



1615 NE EASTGATE BLVD, SEC H
PULLMAN, WA 99163 USA
Phone: 509-592-8144
www.abt360llc.com
alyssa@abt360llc.com
facebook.com/ABT360
instagram.com/ABT_360