

# ABT HOLDING



**ABT HOLDING** medium was developed to provide embryos a sustainable and healthy environment after collection and before transfer or freezing. While ABT Holding is not meant to act as an extended culture medium for incubation, it does contain a variety of components such as amino acids, nutrients and energy sources that facilitate embryo viability.

ABT's Holding starts with ultra-purified Type 1 water. All ABT products are manufactured in a cleanroom environment and go through strict quality control practices, ensuring consistency and quality from lot to lot.

- Includes amino acids, energy sources, Gentamycin as an antibiotic and is HEPES buffered
- 100% traceable components including certified New Zealand sourced BSA
- Filtered to 0.1uM; no additional filtering required
- Tight quality control standards equal consistent results
- Durable EVA packaging with a non-latex needle septum port and removable cap
- Made in the USA with two decades of experience

<u>PRODUCT #</u>	<u>SIZE</u>
<b>ABT024</b>	<b>50 ml Bag</b>
<b>ABT224</b>	<b>20 ml Bag</b>
<b>ABT724</b>	<b>7 ml, Pack of 5</b>

**Storage:** 2 - 8° C  
**Shelf Life:** 18 months from date of manufacture

## ABOUT ABT 360

Located in Pullman, Washington, USA, ABT 360 is owned and operated by Kathy Shearer who has 25 years of experience manufacturing embryo transfer media and Alyssa Hood who has assisted in and managed that same manufacturing facility since 1995. Formed with integrity and experience, ABT 360 has what it takes to support all of your embryo transfer media needs. We look forward to continuing our relationships with the veterinary reproduction community and we will strive to provide quality and consistency in our products, supported by exceptional customer service.



1615 NE EASTGATE BLVD, SEC H  
PULLMAN, WA 99163 USA  
Phone: 509-592-8144  
[www.abt360llc.com](http://www.abt360llc.com)  
[alyssa@abt360llc.com](mailto:alyssa@abt360llc.com)  
[facebook.com/ABT360](https://facebook.com/ABT360)  
[instagram.com/ABT\\_360](https://instagram.com/ABT_360)