

ABT 1 STEP THAW



ABT 1 STEP THAW medium contains 1.0M Sucrose and is intended to be used for the controlled re-hydration of embryos that have been frozen in 10% Glycerol in a single step.

ABT 1 Step Thaw starts with ultra-purified Type 1 water and includes Gentamycin as an antibiotic and certified New Zealand sourced BSA. All ABT products are manufactured in a cleanroom environment and go thru strict quality control practices, ensuring consistency and quality from lot to lot.

- Includes 1.0M Sucrose
- 100% traceable components including certified New Zealand sourced BSA
- Filtered to 0.1uM; no additional filtering required
- Tight quality control standards equal consistent results
- Available in durable EVA bags and 7ml LDPE Single Use Dropper Bottle
- Made in the USA with two decades of experience

<u>PRODUCT #</u>	<u>SIZE</u>
ABT047	50 ml Bag
ABT247	20 ml Bag
ABT747	7 ml, Pack of 5

Storage: 2 - 8° C
Shelf Life: 18 months from date of manufacture

ABOUT ABT 360

Located in Pullman, Washington, USA, ABT 360 is owned and operated by Kathy Shearer who has 25 years of experience manufacturing embryo transfer media and Alyssa Hood who has assisted in and managed that same manufacturing facility since 1995. Formed with integrity and experience, ABT 360 has what it takes to support all of your embryo transfer media needs. We look forward to continuing our relationships with the veterinary reproduction community and we will strive to provide quality and consistency in our products, supported by exceptional customer service.



1615 NE EASTGATE BLVD, SEC H
PULLMAN, WA 99163 USA
Phone: 509-592-8144
www.abt360llc.com
alyssa@abt360llc.com
facebook.com/ABT360
instagram.com/ABT_360